

## Product information


 QF 24 V4  
 V26 2020

# EZ-10 Spin Column COVID-19 Viral RNA Extraction Kit

**Catalog #:** VT 9200  
**Size:** 50 preps, 250 preps, 500 preps  
**Storage:** 4°C to 25°C\*

\*: Product will be shipped at ambient temperature. Once reconstituted, carrier RNA to be stored at -20°C, remaining kit components to be stored at 4°C to 25°C. Components have a one year expiration from date of purchase.

## Product Description:

The kit is **especially designed for COVID-19 (SARS-CoV-2)**. It simplifies isolation of COVID-19 viral RNA from cell-free body fluids with fast spin-column format. No phenol/chloroform extraction is required. Viral RNA binds specifically to the silica membrane while contaminants are removed in the flow-through. PCR inhibitors such as divalent cations and proteins are completely removed in two efficient wash steps, leaving pure viral RNA to be eluted in RNase-free Water. Purified RNA is ready to use in RT-PCR, Northern blotting or other downstream applications.

## Features:

- **Fast:** Using a rapid spin column format, the entire procedure takes about 20 minutes.
- **High Yield:** The recovery yield of viral RNA is generally >85%.
- **Versatile:** Suitable for purification of viral RNA from Nasopharyngeal swabs (dry) or Transport medium/Viral transport medium (VTM) containing Nasopharyngeal swabs, Nasopharyngeal aspirates, and bronchoalveolar lavage (BAL) samples.
- **Non-toxic:** No phenol/chloroform are used.

## Materials Supplied by User:

- Microcentrifuge capable of at least 12,000 × g
- RNase-Free pipettes and pipette tips
- Vortexer
- RNase-Free Ethanol (96-100%)
- RNase-Free Microcentrifuge tubes (1.5 ml or 2 ml)

## Before Starting:

i Freshly Prepare the following **Lysis Mastermix**:

Component	Volume/rxn	Volume/50 rxns	Volume/250 rxns	Volume/500 rxns
Carrier RNA	2 µl	100 µl	500 µl	500 µl x2
Buffer Rlysis-VG	600 µl	30 ml	150 ml	150 ml x2

ii. Mix thoroughly by vortexing.

## Composition:

Components	50 Prep Kit	250 Prep Kit	500 Prep Kit
Buffer Rlysis-VG	30 ml	150 ml	150 ml x2
Universal RPE Solution **	12 ml	30 ml x2	40 ml x3
Nuclease-Free Water***	5 ml	25 ml	50 ml
EZ-10 Spin Columns	50	250	500
2 ml Collection Tubes	50	250	500
Carrier RNA***	1 mg	1 mg	1 mg
Protocol	1		

\*\* : Universal RPE Solution is supplied in a concentrated form. Before use, add 48 ml 96-100% ethanol to 12ml concentrated universal RPE solution and mix well.

\*\*\* : Add Nuclease-Free Water to the tube containing lyophilized Carrier RNA to obtain a solution of 1 µg/µL. Dissolve the carrier RNA thoroughly, divide it into conveniently sized aliquots, and store it at -20°C. Do not freeze-thaw the aliquots of Carrier RNA more than 3 times.

## Storage:

Upon reconstitution, Carrier RNA should be divided into conveniently sized aliquots and stored

at -20°C. Do not freeze-thaw the aliquots of Carrier RNA more than 3 times. All remaining kit components to be stored at 4°C to 25°C. The kit is valid for 1 year.



## Procedure:

### 1. Sample prepara on:

#### A) Nasopharyngeal swabs (dry):

Select between method **a.1** (without Proteinase K)  
OR **a.2** (with Proteinase K - op onal):

**a.1.** Place the swab into a clean 1.5 ml microtube, and snap off the handle. Add 0.2 ml physiological saline, vortex for 30 sec. Then transfer 0.2 ml solu on to a new 1.5 ml microtube and proceed to step 2.

**a.2.** If using Proteinase K (recommended for be er yield):

- i) Add 300 µl of PBS (user supplied) and 10ul of 20mg/ml Proteinase K solu on (user supplied, *Bio Basic cat. #401*) to each swab sample.
- ii) Incubate at 56°C for 5 min. with occasional mixing.
- iii) Centrifuge at 10,000 X g (or max. speed) for 30 sec.
- iv) Use 200 µl from this step as a sample and proceed to step 2.

#### B) Transport medium/Viral transport medium (VTM) containing Nasopharyngeal swabs, Nasopharyngeal aspirates and bronchoalveolar lavage (BAL) samples:\*\*\*

- i. Vortex the tubes containing the swab at max. speed for 1 min.
- ii. Use 200 µl from this step as a sample and proceed with step 2.

\*\*\*Note: If transfer medium is very dilute i.e. if the volume is more than 400 µl, perform the below steps for virus enrichment instead of "i" and "ii":

- iii. Transfer appropriate liquid sample to a new 1.5 ml microtube,
- iv. Centrifuge at 24,000 g for 60 min. at 4°C. Then keep approx. 0.2 ml solu on in the tube but discard the remaining and con nue to step 2.

**2.** Add 0.6 ml of the recently prepared **Lysis Mastermix** into the tube from step 1, vortex vigorously for 30 sec.; incubate at room temp. for 10 min.

**NOTE:** Lysis-Buffer-VG may form precipitate at 4°C, please dissolve it at 65°C and mix well before use.

**3.** Add equal volume of ethanol, mix by inver ng the tube.

**4.** Transfer the mixture into the spin column; keep at room temp. for 2 min.

**5.** Spin at 10,000 g for 1 min., discard the flow-through.

**6.** Add 0.5 ml of Universal RPE Solu on\* to the column, spin at 10,000 g for 1 min. and discard the flow-through.

\*: Universal RPE Solu on is supplied in a concentrated form. Before use, add 48 ml 96-100% ethanol to 12ml concentrated universal RPE solu on and mix well.













**7.** Repeat the Step 6 once.

**8.** Centrifuge at 10,000 g for 1 min., discard the flow-through residue.

**9.** Transfer the column to a new 1.5 ml RNase-free microtube. Add 35-60 µl of Nuclease-Free Water onto the centre of the column; keep at room temp. for 2 min.

**10.** Spin at 10,000 g for 1 minute. Purified viral RNA is ready for use or keep at -20°C for long term storage.

## Index of Symbols

	In Vitro Diagnostic Use		See Instructions for Use
	Catalog Number		Manufacturing Date
	Expiry Date		Number of Tests
	Manufacturer		CE Marking
	Caution		Store at 4- 25°C
	Lot Number		European Authorized Representative



Bio Basic Canada Inc.  
20 Konrad Crescent, Markham,  
ON L3R 8T4, Canada  
Tel: 1-905-474-4493  
E-mail: order@biobasic.com  
Web: www.biobasic.com



Advena Ltd  
Tower Business Centre, 2nd Flr,  
Tower Street, Swatar, BKR 4013, Malta  
Tel: +44 (0) 1926 800 153  
E-mail: info@advenamedical.com  
Web: www.advenamedical.com

